

Public Health Grant

- Ultra cold vaccine freezer \$7,840.00



- New vaccine fridge and freezer \$8,945.33



SO-LOW UNDERCOUNTER FREEZER MODEL MV85-2

**2 Cubic Ft. (54 liters) Upright – Temperature Range: -40°C to -85°C
Solid Door**

Specifications:

* Capacity	2 cubic ft. / 54 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior w x d x h	23" x 27" x 32" 60 x 70 x 81 (cm)
* Interior w x d x h	14" x 18" x 12.5" 36 x 46 x 32 (cm)

PRICE: \$7,840.00
Model MV85-2

Prices above do not include freight
Shipment: 1-2 weeks a/r/o


Features:

Temperature Control	Digital – Battery Backup Datalogger – Hi/Low Alarm
Clearance	4" on sides and back
Doors	Single Solid, locking
Legs	Four, leveling
Insulation	Polyurethane CFC Free
Evaporator	Cold wall
Defrost	Manual
Voltage: Domestic	115/60/1 – 6.3 amps
Supply Plug	NEMA 5-15P
Shipping Weight	265 lbs.
Warranty	Two year parts and labor

Inventory Racks:

Rack # 17-15-2
\$ 115.00 each

Rack holds fifteen
2" high x 5.25" x 5.25"
boxes with dividers.
MV85-2 holds 2 racks.



Specifications subject to change without notice



Quotation 004102761

Follett Products, LLC
801 Church Lane
Easton, PA 18040

Phone: 800-523-9361
 610-252-7301
 Fax : 610-250-0696

Quote-To: 012053
 DAVIS CNTY HOSP
 509 N MADISON ST
 BLOOMFIELD, IA 52537
 US

Bill-To: 012053
 DAVIS CNTY HOSP
 509 N MADISON ST
 BLOOMFIELD, IA 52537
 US

Issue Date
12/16/22 Currency
US Dollars

Project

Spec by:

For Internal Use - 0

Attn: TREVOR COUCHMAN tcouchman@daviscountyhospital.org

Phone: 641-664-2145 EXT 4402

Fax: 641-664-1895

Line	Item Description/Details	Qty	Unit Price	Total
0001	<p>REFFZR16-SS Combination 11.5 cu ft solid door refrigerator; Freezer 4.8 cu ft solid door 16.3 cu ft storage capacity Refrigerator automatic defrost; freezer manual defrost 115V</p> <p>(1) skid @ (337.0 lbs) (34.0 x 29.0 x 79.0 in) Item Reference No.01398536</p>	1	6,686.67	6,686.67
0002	<p>REFFZR21-GS Combination 14.7 cu ft refrigerator, 6.2 cu ft freezer, Low GWP Refrigerant 20.9 cu ft storage capacity Refrigerator automatic defrost; freezer manual defrost 115V</p> <p>(1) skid @ (439.0 lbs) (38.0 x 34.0 x 82.0 in) Item Reference No.01398544</p>	1	8,945.33	8,945.33
***	<p>Please note this quotation number on your purchase order. ***</p> <p>Any freight received damaged must be refused or signed for as damage. Freight claims, repairs, and reimbursement cannot be made after 5 days so Follett Customer Service must be notified of any claim within 48 hours of delivery - See Freight Terms below for details.</p>			

SUBTOTAL 15,632.00

TAX 0.00

US Dollars **Continued**

FCA Point - **FOB Origin, prepay and add**

Cc: MEDTEK RESOURCES 800-397-0222

Please direct all questions to
[Follett Customer Service](#)

Contract Pricing expires on **06/30/23**
 Non Contract Pricing expires on **12/31/22**
 After that date, please contact factory

Please email Equipment POs to
customerservice@follettice.com
 Or fax to 610-250-0696

Unless otherwise specified, freight costs are additional. Lift-gate delivery, at additional cost, will be required if no loading dock is available. Please advise Follett at time of order if lift-gate delivery is needed.



Quotation 004102761

Follett Products, LLC
801 Church Lane
Easton, PA 18040

Phone: 800-523-9361
 610-252-7301
 Fax : 610-250-0696

Quote-To: 012053 DAVIS CNTY HOSP 509 N MADISON ST BLOOMFIELD, IA 52537 US	Bill-To: 012053 DAVIS CNTY HOSP 509 N MADISON ST BLOOMFIELD, IA 52537 US	Issue Date 12/16/22 Currency US Dollars Project Spec by: For Internal Use - 0
Attn: TREVOR COUCHMAN tcouchman@daviscountyhospital.org Phone: 641-664-2145 EXT 4402 Fax: 641-664-1895		

Line	Item Description/Details	Qty	Unit Price	Total
	Also, please see Freight Terms for regulatory limitations of Follett's involvement in damage claims if you contract for your own freight delivery.			

	SUBTOTAL 15,632.00 TAX 0.00 US Dollars 15,632.00
--	--

FCA Point - **FOB Origin, prepay and add**
 Cc: MEDTEK RESOURCES 800-397-0222

Contract Pricing expires on **06/30/23**
 Non Contract Pricing expires on **12/31/22**
 After that date, please contact factory

Please direct all questions to
[Follett Customer Service](#)

Please email Equipment POs to
customerservice@follettice.com
 Or fax to 610-250-0696

Unless otherwise specified, freight costs are additional. Lift-gate delivery, at additional cost, will be required if no loading dock is available. Please advise Follett at time of order if lift-gate delivery is needed.

This quote contains confidential pricing information intended only for the use of the individual or entity named above and may be legally privileged. If you are not an employee or agent of the intended recipient, you are hereby notified that any disclosure, copying, or distribution, in part or in whole, to any party not associated with the acceptance of this specific quote, or the taking of any action based on the contents of this quote, is strictly prohibited.

Warranty Summaries

Unless otherwise specified, all warranties cover repair or replacement of the product by an authorized service company.

Symphony Plus ice and water dispensers and Maestro Plus ice makers

- 3 year parts and labor on entire machine
- 5 year parts on compressor

Horizon ice machines/ Ice Manager Diverter Valve System

- 3 year parts and labor on entire machine
- 5 year parts on compressor

Slope front bins, carts and accessories

- 3 year parts
- 1 year labor

Vision undercounter dispensers

- 2 year parts and labor

Ice Pro semi-automatic bagging systems

- 2 year parts and labor

Upright bins, Ice-DevIce, and ITS storage and transport systems

- 5 year parts and labor

7 Series / 15 Series Dispensers

- 2 year parts and labor on entire machine
- 5 year parts on compressor
- For indoor use only. Designed for commercial use. Follett is not able to provide in-house service for residential installations

Countertop and undercounter refrigerators and freezers

- 2 year parts and labor on entire machine
- 5 year parts on compressor

Upright refrigerators and freezers

- 2 year parts and labor on entire machine
- 5 year parts on compressor

Warranty Statements

- For full warranty terms, conditions, and exclusions please refer to the Warranty Statements available at <http://follettice.com/support/>.

For warranty service, contact Follett's Technical Service department by telephone or email –

- 877.612.5086 – US, Canada, Mexico
- 610.252.7301
- Techsupport@follettice.com

Terms

- This price list supersedes all others
- Net 30 day terms on approved accounts
- Prices subject to change without notice
- Specifications subject to change without notice

Freight Terms and Freight Damage

Freight is an additional charge. Standard delivery rates include dock-to-dock services. Special services such as, but not limited to, lift gate or inside delivery are not included. Please advise Follett LLC at time of order if special services are needed.

All freight is to be inspected upon delivery. If visible signs of **damage** exist, please refuse delivery or sign your delivery receipt “damaged.” **Follett Customer Service must be notified within 48 hours.** Wherever possible, please include detailed photos of the damage with the original packaging so that we may start the freight claim process.

Freight industry policies prohibit Follett from taking responsibility for freight claims where Follett has not contracted with the freight carrier. Follett recommends shipment with one of our contracted carriers and believes we can provide these services at a reasonable cost to you. If you choose to contract for the delivery of your equipment using your own 3rd party carrier (even if it is a carrier that Follett uses), Follett will be unable to submit claims for freight damage or provide restitution for freight damage expenses. Should you experience freight damage with your carrier, Follett will work with you to document the event so that you can file claims with your carrier.

Please contact Follett Customer Service for assistance with ordering replacement parts and/or equipment. Equipment and parts damaged due to improper handling and/or receiving may void the Follett warranty provision. Any damaged (cosmetic or otherwise) parts/equipment placed into service must be reported to and approved by Follett to honor warranty claims.

Return Policy

Follett equipment may be returned under the following conditions:

1. The equipment is new, unused, and not installed.
2. A return materials authorization (RMA) accompanies the returned equipment. Such RMA must be requested from Follett Customer Service within 45 days after shipment
3. Equipment must be received at the Follett factory in Easton, PA within 30 days of the issuance of the RMA by Follett
4. There is a 20% re-stocking fee for all returns after equipment is inspected at Follett. If the cost to return the equipment to stock exceeds 20%, the actual cost will be deducted
5. Return freight is the responsibility of the customer returning the equipment.
6. Units to be returned must be in brand new condition and in the original cartons or boxes. If you do not have the original packaging, please contact Follett Customer Service to order new packaging. If returned equipment is damaged because of improper packaging, Follett will not be held responsible.

Please Note: Special quote request (SQR) items are not returnable to Follett.

Help us help you keep your order on track

To avoid delays, please include the following information on your purchase order:

- Quotation number
- Purchase order number
- Billing address
- Shipping address
- Freight terms (if different than quotation)
- Special delivery instructions (required if not standard dock-to-dock delivery)
 - Lift gate
 - Inside delivery
 - Follett authorized delivery or installation service
- Pricing that matches quotation
- Model numbers that match quotation
- Net 30 payment terms (for approved accounts)
- Tax exempt/resale certificate

Please contact Follett Customer Service for assistance

Hours Monday – Friday 8:00 am to 5:00 pm (Eastern Standard Time)
 Phone 1.610.252.7301
 Toll free 1.800.523.9361

Click here to view installation and service videos that will reduce your on-site time and ensure a successful installation/service call.

MPR-715F-PA | REFFZR21-GS



MPR-715F-PA | REFFZR21-GS

Pharmaceutical

Refrigerator with Freezer
EPA SNAP Compliant Refrigerants

+2°C to +14°C Refrigerator
-20°C to -30°C Freezer

14.7 cu.ft. | 415 L
6.2 cu.ft. | 176 L

This space-saving cabinet combines both refrigerator and freezer functions within a single footprint. It is designed for use as a pharmaceutical or vaccine refrigerator or a biomedical or general use freezer. The slim profile permits quick access to all products. Because less lab floor space is required, this model can be located in narrow hallways.

Independent Cooling Systems

Separate refrigeration systems create independently controlled chambers. The freezer system maintains temperatures as low as -30°C. The refrigerator compartment operates at temperature setpoints over a range of +2°C to 14°C. The integrated refrigeration system is engineered for quiet operation.

Defrost

The refrigerator includes a unique automatic defrost system which activates electronically when needed. The refrigerator operates above freezing at all times. This prevents sensitive vaccines, reagents and other materials from accidental freezing. The freezer can be defrosted manually when required.

Superior Uniformity

Forced air circulation through a back wall plenum in the refrigerator provides even air distribution for precise temperature uniformity and efficient temperature recovery to setpoint quickly following routine door openings. The freezer is configured with a cold wall cooling system with external coils strategically positioned within the insulation for excellent uniformity.



Sustainable Refrigerants

Each of the two refrigeration circuits uses environmentally friendly refrigerants that will not compromise reliability or operational performance. These refrigerants enable the system to recovery setpoint temperature quickly following door openings, and to resist high ambient temperatures.



Microprocessor Control Panel

The microprocessor temperature controller maintains accurate temperature at desired setpoint, and responds quickly to changes in ambient temperature and door openings. Audible and flashing LED visual alarms alert the user to any temperature deviations from setpoint. Both refrigerator and freezer temperatures are independently displayed on the control panel.



Access, Storage and Security

The refrigerator includes two insulated doors with triple pane glass viewing windows. Heat reflection film minimizes condensation from ambient humidity. The freezer has insulated solid doors. A keyed door lock assures inventory security. Refrigerator has three adjustable wire shelves and the freezer has two wire shelves.

Pharmaceutical Refrigerator with Freezer

MPR-715F-PA | REFFZR21-GS

SNAP Refrigerants and Compliance with EPA Recommendations


This model operates on new generation SNAP compliant, low flammability HFO refrigerants to assure low environmental impact. The use of SNAP compliant refrigerants (Significant New Alternatives Policy) established by the United States Environmental Protection Agency affirms our commitment to environmental sustainability by integrating the latest refrigeration technology. For information on SNAP refrigerants visit <https://www.epa.gov/snap>

Validated Storage

The PHCbi pharmaceutical refrigerator freezer combination is designed to assure stored product safety to meet CDC and third-party guidelines for vaccine and biological material storage.

Refrigerator Defrost Temperature Data

The automatic defrost system eliminates temperature spikes to maintain safe storage conditions. Stored products see stable temperature at all times.



Dimensions

Unit : inch | mm




Illustration shown for Model MPR-715F-PA | REFFZR21-GS. Dimensions apply to models with windows and models with solid doors.

Model Number		MPR-715F-PA REFFZR21-GS refrigerator / freezer
External Dimensions (W x D x H) nominal	inches mm	35.4 x 28.1 x 75.2 900 x 715 x 1910
Internal Dimensions (W x D x H) nominal	inches mm	31.9 x 24.2 x 35.2 810 x 615 x 894 15.1 x 21.7 x 16.6 385 x 552 x 422 (each chamber)
Volume	cu.ft. liters	14.7 415 / 6.2 176
Net Weight	lbs. kg	357 162
Chambers	qty	1- Top / 2- Lower
Performance		
Warranty ¹⁾		2 years parts and labor, 3 years compressor parts
Temperature Control Range	°C	+2 to +14 / -15 to -30
Vaccine Storage Operating Temperature	°C	+2 to +8 / -15 or colder
Factory Pre-Set Temperature	°C	+5 / -30
Highest Ambient Temperature and Maintains Cabinet Temperature	°C	+30
Protects Vaccine from Freezing		Operates above freezing at all times
Control		
Microprocessor Controller, Adjustable	°C	Increments of 1 - Independently controls refrigerator & freezer
Digital Temperature Display		LED
Controller Security		Lockable with key pad
Electronics Diagnostics		Total control system
Refrigeration		
Cooling Method Internal Airflow for Precise Temperature		Forced air / Direct cold wall cooling
Defrost Method		Automatic ²⁾ / Manual
Refrigeration System	qty	2- Independent Air cooled, CFC free
Refrigerant ³⁾		SNAP Approved HFO R-513A
Insulation		CFC free urethane
Construction		
Outer Door, Swing door with CFC Insulation	qty	2-Top, Bi-parting with dual pane glass 2-Bi-parting, Solid
Interior		Painted steel / Painted steel
Exterior		Zinc galvanized steel, Acrylic finish
Outer Door Lock	qty	1-Key, Center left, Locks top left and lower left 1-Key, Center right, Locks top right and lower right
Interior Light	qty	1- LED, control panel switch / 0
Shelves	qty	(3) Adjustable, Wire / (2) Wire
Casters	qty	2-Swivel, 2-Fixed
Adjustable Feet	qty	2- Front of base, for securing unit in place
Access Port	qty	2- Side wall, Left top and bottom / 0
Access Port Diameter	inches mm	1.1 30
Alarms (V = Visual, B = Buzzer, R = Remote)		
Power Failure		V-B-R
High Temperature		V-B-R
Low Temperature		V-B-R
Door Open		V-B
Remote Alarm Contacts		Normally open, Normally closed, Common
Remote Alarm Output		DC 24V 2A
Electrical and Noise Level		
Power Supply		115V, 1Ø, 60Hz, NEMA 5-15P requires 5-15R receptacle
Noise Level ⁴⁾	dB(A)	43

¹⁾ Current warranty offered at time of printing and may be subject to change; US and Canada only

²⁾ Electronically monitored defrost only when needed. Protects vaccines from freezing

³⁾ Compliant with SNAP Acceptability Determination https://www.epa.gov/sites/production/files/2016-09/documents/snap_notice_31_factsheet.pdf

⁴⁾ Nominal value, background noise 20 dB(A)

Specifications are subject to change without notice.

Performance data herein is based on independent testing at time of publication.



PHC Corporation of North America